

# Widely tunable CW Optical Parametric Oscillators: Mastering the challenges posed in quantum technology research

## Korbinian Hens<sup>1)</sup>, Jaroslaw Sperling<sup>1)</sup>, Niklas Waasem<sup>1)</sup>, Ronja Gärtner<sup>1)</sup>, Gunnar Elgcrona<sup>2)</sup>

### Introduction

Widely tunable continuous wave optical parametric oscillators (cw OPOs) are gaining popularity as novel sources of tunable laser light, not least due to the unprecedented wavelength coverage in the visible and the near infrared spectral range. While the technology is becoming increasingly recognized to offer a number of potential advantages in particular to the quantum research experimental community, the requirements are often challenging the device performance. In this context, we discuss the characteristics of state-ofthe-art tunable cw OPO designs and describe several tuning schemes tailored to meet various experimental needs.

In an illustrative fashion, we compare recently published experimental data photoluminescence several excitation (PLE) studies, which have been carried out on ensembles as well as on individual quantum emitters.

### **OPO Characteristics**

The huge potential of OPOs, compared to conventional sources of laser light, derives from their exceptional wavelength versatility, as they are in principle not limited by the wavelength coverage of suitable laser gain media<sup>1</sup>. Only relatively recently, advances in the design of nonlinear crystals and the availability of high performance solid-state pump laser sources enabled the practical realization of cw OPOs with game-changing characteristics.

### **C-WAVE** | The tunable laser light source

- Broad tuning range
- IR (900 1300 nm)
- VIS (450 650 nm)
- High output power
- IR: 400 mW, typ. max. > 1.2 W
- VIS: 200 mW, typ. max. > 600 mW
- Single-frequency operation
- Linewidth < 1 MHz (typ. < 500 kHz)</p>
- Fast mode-hop-free tuning: > 20 GHz , > 15 GHz / s (VIS)
- All solid state no consumables

### **PHOTOLUMINESCENCE EXCITATON STUDIES OF SINGLE-PHOTON EMITTERS**



Photoluminescence excitation spectra of an ensemble of Si-V centers recorded at two excitation powers (a) and of an ensemble of Ge-V centers (b) are recorded at room temperature. The data has been recorded by exploiting the full coarse tuning range of the C-WAVE OPO with a step-size of approx. 10 nm, straight-forward accessible by automated crystal selection and temperature tuning. Adapted with permission from reference 2.



A further refined electronic level structure of Ge-V can be investigated at cryogenic temperature conditions at around 602 nm.



PLE spectra of transition C (left) and D (right) are recorded under tuning resonant excitation using the truly continuous (mode-hop-free) tuning C-WAVE. Zero detuning mechanism of corresponds to 602.2903 nm and 602.4828 nm respectively. Orange lines represent the signal when an additional 532 nm gating laser is switched on, blue lines represent the signal when the gating laser is off. Adapted with permission from reference 3.

<sup>1</sup> Division Hübner Photonics, HÜBNER GmbH & Co. KG, Heinrich-Hertz-Str. 2, 34123 Kassel, Germany <sup>2</sup> Division Hübner Photonics, Cobolt AB, Vretenvägen 13, 17154 Stockholm, Sweden





The OPO resonator modes (grey color) undergo a coarse pre-selection by the thin Etalon modes (magenta color). In combination with the thick Etalon mode (blue color) only a single OPO resonator mode (red color) is supported for resonant oscillation contributing to the laser light output.

laser

2 2.1 2.2 2.3

Energy E (in eV)

1.8 1.9

<sub>L</sub> = 2.29 eV

 $E_{\perp} = 2.09 \text{ eV}$ 



**Excitation Detuning (in GHz)** 



A very recent study<sup>4</sup> on defect centers in hexagonal boron nitride (hBN) uses the tuning capabilities of C-WAVE for characterization at room temperature. The figure shows a photoluminescence spectrum of a single defect center in hBN.

Blue shaded area: integration interval over zero phonon line (ZPL) for PLE spectrum, yellow shaded area: integration interval over longitudinal optical (LO) sideband for PLE spectrum,

blue line and yellow line: corresponding PLE spectra.

Photoluminescence spectra of the single defect center in hBN are recorded under laser excitation energy tuning from 2.09 eV (593.30 nm) to 2.29 eV (541.48 nm) by employing a suitable crystal temperature tuning protocol of the cw OPO that is used for the emitter excitation. The shown data is corrected for background and Raman-lines. Adapted with permission from reference 4.

	323.710
	323.700
	323.690
ŦŦ	323.680
cy /	323.670
luen	323.660
Free	323.650
	323.640
	323.630
	323.620

the Thick Etalon.



wavelength range.

### Summary

- OPO technology provides an unprecedented wavelength coverage
- The tuning capabilities of C-WAVE cover the relevant frequency scales for quantum technology research
- The concept presented is quite general and thus, it allows for continuous adaption to new experimental requirements respectively novel types of singlephoton emitters and systems alike
- Overall, cw OPO technology is about to mature into a recognized choice among laser light sources that advance the rapidly evolving field of quantum research

### References

- [1] Ebrahim-Zadeh, M., [Optical Parametric Oscillators], Handbook of Optics 2nd Ed. Chap. 22, McGraw-Hill, 22.1-22.72 (2001)
- [2] Häußler, S., Thiering, G., Dietrich, A., Waasem, N., Teraji, T., Isoya, J., Iwasaki, T., Hatano, M., Jelezko, F., Gali, A. and Kubanek, A., "Photoluminescence excitation spectroscopy of SiV- and GeVcolor center in diamond", New J. Phys. 19(6), 063036 (2017).
- [3] Chen, D., Mu, Z., Zhou, Y., Fröch, J. E., Rasmit, A., Diederichs, C., Zheludev, N., Aharonovich, I. and Gao, W., "Optical gating of resonance fluorescence from a single germanium vacancy color center in diamond", Phys. Rev. Lett. 123(3), 033602 (2019)
- [4] Wigger, D., Schmidt, R., Del Pozo-Zamudio, O., Preuß, J. A., Tonndorf, P., Scheider, R., Steeger, P., Kern, J., Khodei, Y., Sperling, J., de Vasconcellos, S. M., Bratschitsch, R. and Kuhn, T., "Phononassisted emission and absorption of individual color centers In hexagonal boron nitride", 2D Mat. 6(3), 035006 (2019)

### Acknowledgement

The authors gratefully acknowledge the support and fruitful discussions with Alexander Kubanek and his group, Wei-bo Gao and his group, and Rudolf Bratschitsch and his group.

Corresponding author: Dr. Korbinian Hens, Product Manager – Korbinian.Hens@hubner-germany.com



### **HÜBNER Photonics**