

# Novel narrow linewidth 785 nm diode laser with enhanced spectral purity facilitates low-frequency Raman spectroscopy


Magnus Rådmark, Gunnar Elgcrona, Håkan Karlsson,  
Cobolt AB, a part of HÜBNER Photonics, Solna, Sweden

## INTRODUCTION

Raman Spectroscopy enables fast, sensitive and label-free chemical analysis of a large range of materials and has become a routine analytical tool in a wide range of material science and process-control applications.

As the Raman signal is weak it is critical that the illumination laser has a very high level of spectral purity, for efficient detection of the Raman signal. Most materials can be characterized by studying Raman shifts down to  $100\text{ cm}^{-1}$ , but in some cases, for instance for determining the crystallinity of pharmaceutical compounds, it is required to study Raman shifts in the low-frequency regime;  $<100\text{ cm}^{-1}$ .

Recording Raman signals in the low-frequency regions is possible using VBG based notch filters, that can provide blocking of  $>60\text{ dB}$  at as close as  $5\text{ cm}^{-1}$  to the laser peak, and by using laser sources with a very high level of spectral purity and stability.




## LASER REQUIREMENTS

In addition to a stable narrow linewidth and low wavelength drift over both time and ambient temperature, which is required for high spectral resolution, the laser used for low-frequency Raman needs to have a high level of spectral purity (or low Side-Mode Suppression Ratio, SMSR) very close to the main laser peak.

Laser requirements for Low-frequency Raman:	
Linewidth	$<\sim 30\text{ pm}$
Wavelength stability	$<10\text{ pm}$ (over min)
Spectral purity	$>60\text{dB}$ @ $<5\text{cm}^{-1}$ ( $\sim 300\text{ pm}$ )

An illumination wavelength around  $785\text{ nm}$  is preferred for many materials due to an optimum compromise between Raman signal strength, auto-fluorescence level and detector sensitivity.




Conventional wavelength-locked diode lasers need external filtering with an additional VBG clean-up filter in order to meet the requirements on spectral purity for low-frequency Raman due to broad-band Amplified Spontaneous Emission (ASE) from the diode.

## TECHNICAL SOLUTION

Patent pending

To overcome the challenge with broad-band ASE from wavelength-locked diode lasers we developed a design where we use a highly reflective VBG for locking instead of a partially transmissive one, as in the conventional case.

An intra-cavity polarizing element and a polarizing beam splitter are used to control the level of feed-back from the VBG to the emitter and the output coupling out of the cavity. The output-coupled emission is already cleaned-up from ASE as the ASE still leaks out from the locking VBG. In this way, over  $60\text{ dB}$  spectral purity at  $<300\text{ pm}$  from the main peak is achieved without any further external filtering.




$>60\text{dB}$  spectral purity @  $<300\text{ pm}$  without external filtering



To achieve stable locking and stable spectral performance over time and temperature we used HTCure™ technology for the laser manufacturing; all optical elements are mounted onto a single temperature-stabilized platform and cured with high-temperature adhesives. The optical cavity platform is then mounted into a hermetically sealed package.





## PERFORMANCE




SMSR of  $>64\text{ dB}$  at  $\pm 300\text{ pm}$  (corresponding to  $\pm 5\text{ cm}^{-1}$ ) and  $>74\text{ dB}$  at  $\pm 1\text{ nm}$  over  $15\text{-}55\text{ °C}$  temperature cycling



Linewidth (FWHM) of  $<31\text{ pm}$  (corresponding to  $0.5\text{ cm}^{-1}$ ) in temperature cycling over  $15\text{-}55\text{ °C}$  and  $<27\text{ pm}$  at constant temperature.



Wavelength drift  $<4\text{ pm}$  at constant temperature and  $<17\text{ pm}$  in temperature cycling over  $15\text{-}55\text{ °C}$



## REFERENCES

- M.P.T. Fraser, et al. "Use of low-frequency Raman spectroscopy and chemometrics for the quantification of crystallinity in amorphous griseofulvin tablets" Vibrational Spectroscopy 2015, 77, 10-16.
- P.J. Larkin, et al. "Polymorph Characterization of Active Pharmaceutical Ingredients (APIs) Using Low-Frequency Raman Spectroscopy". Applied Spectroscopy. 2014, 68 (7):758-776
- A. Glebov, et al. "Volume Bragg Gratings as Ultra-Narrow and Multiband Optical Filters" Proc. of SPIE Vol. 8428 84280C-1