## Laser and **Terahertz Technology**

HÜBNER Photonics | Coherence matters.

**HÜBNER Photonics** 



#### HÜBNER Group

Mobility. Materials. Photonics. | united by pass

#### HÜBNER Photonics, a part of the HÜBNER

Coherence matters.

 $(\mathbf{i})$ 

Ζ

ш

Ζ

Laser and Terahertz technology

A bright way into the future

#### Laser technology

High performance lasers

- Cobolt Single and multi-line lasers
- C-WAVE Tunable lasers
- C-FLEX Laser combiners
- VALO Femtosecond lasers

Product overview

Applications

#### Terahertz technology

Visualizing the invisible

#### Contact

HÜBNER Photonics worldwide

Mobility. Materials. Photonics. | united by passion.

sion.	4
Group	
	6
	8
	10
	10
	12
	14
	20
	24
	26
	28
	32
	40

44

## HÜBNER Group Mobility. Materials. Photonics. united by passion.

Since the beginning in 1946 when Kurt Hübner founded the company, Today the HÜBNER Group consists of more than 30 different com-HÜBNER has been shaped by a spirit of enterprise, inventiveness and panies and has more than 3,500 employees around the world. And HÜBNER products continue to stand out with their high-quality workinnovation. Today it continues to be this same passion for new solutions that drives the company and its employees. Again and again over manship, their superior convenience and their distinctive durability. the years, new concepts for new products and new markets have been This high level of quality is to be found in all of the Group's business developed and implemented. sectors and throughout the entire product range.

It is on this basis that HÜBNER has become the company that it is This all comes together in our company slogan: today: A comprehensive system supplier for transportation technology and the worldwide leader for gangway systems (Mobility), a well-recognized specialist for sophisticated elastomer solutions and products, insulation and composite materials (Material Solutions), as well as an up-and-coming company in the field of laser and terahertz technologies (Photonics).



Mobility. Materials. Photonics. | united by passion.

## HÜBNER Photonics | Coherence matters.

Coherence matters in our customer relations - we operate in phase with our customers' needs and therefore stand out by maximizing customer value.

"Two wave sources are perfectly coherent if they have a constant phase difference and the same frequency." In a way, this condition sums up everything HÜBNER Photonics stands for.

We not only make game-changing lasers and light sources, but also rethink all kinds of other wave technologies including terahertz imaging and spectroscopy. The proven corporate values of the HÜBNER Group are brought together with innovative ideas and topnotch technologies.

Coherence for us is far more than the perfect interaction of waves. It is the basis of our daily work and the fundament on which we build our business and our relations within our company and with our partners.

Coherence matters for our products – since we are committed to supplying reliable, innovative high performance photonics solutions.

Coherence matters in our attitude - our teams are built up by worldleading expertise in photonics, moving as one.

Coherence matters is in our genes and our spirit – all day, every day.





High-performance lasers for science and industry HÜBNER Photonics offers a wide assortment of high performance laser products which meet the ever increasing opportunities for lasers in science and industry. This portfolio includes compact single frequency CW lasers, single frequency tunable CW lasers, multi-line lasers, diode lasers, Q-switched lasers, femtosecond lasers and laser combiners covering the full UV-Visible-MidIR spectrum. Main application markets are to be found in spectroscopy, holography, quantum optics and many others. The combination of competencies across the HÜBNER Photonics division facilitates volume manufacturing of exceptional laser designs with unparalleled reliability. This means that life science and analytical instrument manufacturers now have a supplier they can really rely on and research scientists have a reputable source for lasers to meet their needs.

## TERAHERTZ Terahertz technology

An innovative tool for the detection, analysis and imaging of materials Outstanding ease of use, extremely compact and entirely safe for persons and materials - these are the characteristics of the terahertz systems developed by HÜBNER Photonics on the basis of cutting edge research. Using terahertz waves, the new devices make it possible to visualize concealed contents in various situations. For example, in letters and packages, suspicious substances and objects can be detected immediately and analyzed and visualized with hyperspectral imaging within seconds. The HÜBNER Photonics imagers and spectrometers are therefore especially well suited for use in mailrooms and other settings where security can be an issue. The devices also offer unique capabilities for non-destructive quality control in connection with industrial production.



Introduction of Cobolt Skyra™ and C-FLEX

#### Acquisition of Cobolt AB

1952

**Diversification into laser technology** in collaboration with Fraunhofer Institute

2008

2006

Diversification from HÜBNER into terahertz technology in collaboration with Fraunhofer Institute

Production of rubber tubing connections as gangways for railway passenger cars of the Deutsch Bahn and production of the folding bellows for articulated buses

1946 Founding of the company by Kurt Hübner in Kassel, Germany. Start of rubber production activities

## History A bright way into the future

Again and again, HÜBNER has been a driving force for technological progress. The family-owned company places a high priority on the innovative value and quality of its products and concepts - from the initial idea to the market-ready product.

In the field of laser and terahertz technologies, HÜBNER made a commitment to basic research early on and began to cooperate with leading scientific institutes. The company secured further know-how through its integration of the Swedish laser manufacturer Cobolt in 2015.

Even before the company's activities had been combined into the new business segment HÜBNER Photonics, two in-house developments received the renowned photonics Prism Award - the tunable laser C-WAVE as well as the terahertz spectrometer T-COGNITION.

In addition to its involvement in basic research, HÜBNER Photonics is also committed to ensuring the utility of its developments for customers in the real world. Many products and solutions are developed in close cooperation with industrial customers and with an eye on their specific needs. And these products also undergo continual enhancement and improvement based on the demands of commercial use. In addition, all customers benefit from the wide cross-industry know-how of the HÜBNER Group.



2000 Founding of Cobolt AB

by Cobolt

Lasers have long played a special role in Kassel, where the HÜBNER headquarters is located. For more than 40 years, green laser lines have marked the night sky of the city. Kassel's "Laserscape" created for the Documenta exhibition in 1977, is the world's first permanent laser light artwork.

# LASER TECHNOLOGY

Cobolt

## **High-performance lasers** A wide array of laser sources

Laser technology has come to have an indispensable role in scientific and medical research, in industrial production and even in many everyday applications from cosmetics to autonomous vehicles. At HÜBNER Photonics, our laser technology is offered in 4 product lines:



**Cobolt – Single and multi-line lasers** Through the well-known Swedish laser manufacturer Cobolt, a proven supplier of high-performance lasers of more than 20 years, HÜBNER Photonics offers one of the industry's broadest ranges of compact single-frequency CW lasers, diode lasers and pulsed lasers across the full UV-Visible-MidIR spectrum.



C-WAVE - Tunable laser Complementing Cobolt's broad portfolio of compact lasers, the tunable lasers developed by HÜBNER Photonics cover an unusually wide spectrum and have gained a strong reputation throughout the industry. C-WAVE is one of the most unique kinds of lasers on the market - a single-frequency, CW, frequency doubled OPO (optical parametric oscillator), providing broadly tunable emission across the visible (and NIR) spectrum at the click of a button.



**C-FLEX – Laser combine** With the trend towards simplicity and user friendliness, the C-FLEX laser combiners offer the flexibility of combining any lasers from the extensive range of high-performance lasers from Cobolt.



VALO – Femtosecond fiber lasers A natural compliment to the high-performance laser portfolio is the VALO ultrashort fiber lasers. Thanks to an innovative laser design with sophisticated control of dispersion and non-linear effects, the VALO lasers offer ultra-short pulses from a compact, robust and turn-key system.





## Cobolt -Single and multi-line lasers

Innovative laser solutions with high output power, stable single mode operation and large wavelength flexibility in the UV-Visible-MidIR spectral range.

Cobolt lasers have become a preferred choice by leading instrument manufacturers and scientists in some of the most challenging applications in the fields of biomedical research, on-line quality and process control, clinical diagnosis, material research, particle analysis and semiconductor metrology. Our commitment is to increase the availability of laser-based equipment that can contribute to improving the quality of life and a better environment.

A combination of sophisticated laser designs and the HTCure<sup>™</sup> Technology for advanced laser manufacturing allows Cobolt to provide lasers in compact formats with the performance and reliability required by the most demanding applications and for use in industrial environments.

## Applications

- Flow cytometry
- Interferometry
- Semiconductor metrology
- Gas detection
- Materials processing
- Optogenetics
  - Optical tweezers
  - Holography
  - LIBS LIDAR

14

- Fluorescence microscopy
- Super-resolution microscopy
- DNA sequencing and analysis
- Raman spectroscopy

Dynamic light scattering

## **Cobolt** CW and pulsed lasers

Compact single-frequency CW lasers, diode lasers and Qswitched lasers across the full UV-visible-MidIR spectrum. Using proprietary HTCure<sup>™</sup> laser manufacturing technology, Cobolt lasers are associated with outstanding reliability and lifetime.



		1105-01-5-	1 Seites	1 Seiles	d'a	per sole consister					
	cobolt Od .	cobolt 06	Copolt 08	cobolt Ski	Copot PC	copoit of	copolitio				
			<b>N</b>					max. power level			
CW lasers	•	٠	•	٠	•			3000 mW			
Pulsed lasers						٠	٠	1000 mW			
Single frequency	•		•					3000 mW			
Modulated lasers		٠		٠				300 mW			
Narrow linewidth			٠					500 mW			
Tunable lasers						•		80 mW			
Multi-line lasers				•				100 mW			
Fiber coupled/pigtailed	•	٠	٠	٠							

#### **CW** lasers

#### Cobolt 04-01 Series and 05-01 Series

#### $\bullet \bullet \bullet$

#### Powerful single frequency CW diode pumped lasers

- 320 1064 nm up to 3 W
- < 0.1 % rms noise, in a perfect TEM00 beam</p>
- < 1 MHz linewidth, superior spectral purity and wavelength stability

#### Cobolt 06-01 Series

#### 

- 375 1064 nm up to 400 mW
- Modulated laser diode (MLD)) and diode pumped lasers (DPL) with fully integrated electronics
- Fast and deep modulation, fiber pigtailed option

#### Cobolt 08-01 Series

#### 

#### Compact narrow-linewidth lasers

- 405 1064 nm up to 500 mW
- SLM diode pumped lasers (DPL) and frequency stabilized diode lasers (NLD)
- Integrated optical isolators/filters

#### Cobolt Skyra™

#### 

- Up to 4 laser lines, 405 660 nm
- Permanently aligned in a single beam
- Fully integrated electronics
- Fiber coupled option



### Cobolt Rogue™ Series

High power, CW diode pumped lasers

- 640 nm, up to 1 W
- <150 GHz spectral bandwidth</p>

#### **Pulsed lasers**

## Cobolt Tor<sup>™</sup> Series

High performance Q-switched lasers

- 355 nm, 532 nm, 1064 nm
- <5 ns, >7 kHz free running
- Single shot to 1 kHz triggerable
- up to 0.5 mJ/pulse

#### Cobolt Odin™ Series Compact, tunable Mid-IR OPOs

- Wavelength selectable 3,1 µm 4,6 µm
- Tunable up to 50 nm, narrow linewidth option
- Up to 80 mW at 10 kHz





## **Outstanding** from start to finish



#### Pioneer with a patent

The Swedish company Cobolt arose from a research project in the year 2000. Today as part of HÜBNER Photonics, it is one of the world's leading developers and producers of high-performance lasers. The product portfolio encompasses innovative laser solutions with high output power and stable single-mode operation. In the UV-Visible-MidIR spectrum range, the products offer exceptional flexibility in wavelength adjustment.

Thanks to the sophisticated cavity designs and patented HTCure<sup>™</sup> production technology, Cobolt can produce compact lasers for a wide range of applications as well as for industrial use. Providing dependable, reliable, long-lasting laser technology for better quality of life and improved environmental protection – that's the aim of our company.

#### The HTCure<sup>™</sup> process

HTCure<sup>™</sup> is a high-temperature curing process developed by Cobolt as part of the laser production process. It allows for a very high degree of reliability in the curing process, exceptional optical performance and unparalleled stability.

In the HTCure<sup>™</sup> process, all of the components of the assembly are selected for their thermomechanical stability and the lasers are packaged in a hermetically sealed sub-housing. As a result, during the production the complete laser unit may be heated to over 100 degrees Celsius for several hours and in various steps without losing its alignment or suffering any damage.

This thermal curing step ensures an extremely rigid and at the same time robust solution so displacements do not occur later on.







## **C-WAVE** The tunable laser light source

C-WAVE is the tunable laser light source for continuous-wave (cw) Focus on your research, not on laser handling: C-WAVE helps you free emission in the visible and near-infrared wavelength range. Its techyour mind for your main tasks. nology is based on optical parametric oscillation (OPO) and it is fully computer controlled. Thus, it allows you to tune from blue to red and **Applications** into the near-infrared without any change of dyes or optical components. This makes C-WAVE a flexible and user-friendly laser for your Quantum optics applications.

#### Change the way you work

Visible, widely tunable, continuous-wave - for a long time this meant handling of toxic dyes, changing laser media or resonator mirrors and often a restricted narrow tuning range. C-WAVE is a solid state system that has no consumable components such as dyes. The wavelength can be simply set at the computer.

C-WAVE tunes itself automatically and guarantees superior beam quality as well as output stability across the whole tuning range offering both high flexibility and precision at the same time. It provides single frequency operation, narrow spectral linewidth and options for frequency stabilization combined with an unprecedented spectral **Raman spectroscopy** coverage.



- Holography
- Nanophotonics
- Atomic physics
- Metrology
- High resolution spectroscopy
- Biophotonics
- Photochemistry
- Interferometry

## **C-WAVE** Flexibility – with precision



#### **Operation principle**

C-WAVE combines two nonlinear processes to achieve its outstand-The concept of C-WAVE allows the user to build laser-light sources for tunable continuous-wave emission from the near UV to the ing spectral coverage: In the first step (OPO), a fixed frequency laser pumps a nonlinear periodically poled crystal. Signal and idler phoinfrared. tons with tunable frequencies in the near-infrared wavelength regime are generated. Output power – tailored to your needs

Subsequent second harmonic generation (SHG) using a frequency doubling crystal leads to conversion of the tunable infrared light into the visible wavelength range.

The modular design of C-WAVE enables you to choose the color ranges that you require!



#### **AbsoluteLambda™**

External wavelength measurement device (wavemeter) and software package to provide enhanced absolute wavelength accuracy as well as precise frequency stabilization for C-WAVE.





Exemplary measurement of stabilized wavelength using AbsoluteLambda.



#### Wide tunability – choose your colors

Depending on the required output power level, C-WAVE is either pumped by an external single-frequency laser or comes with an integrated laser, making operation and application even easier for you.

Typical output power of different C-WAVE models

375 nm 395 nm 405 nm • 415 nm 425 nm 445 nm 457 nm 473 nm 488 nm 491 nm 505 nm 515 nm 520 nm 532 nm 553 nm 561 nm 594 nm • 633 nm 638 nm 640 nm 647 nm 660 nm 685 nm 690 nm 705 nm 730 nm 760 nm 785 nm 808 nm 808 nm 852 nm 915 nm 940 nm 975 nm

1064 nm

## C-FLEX Laser combiners Choose your colors

The C-FLEX laser combiner harnesses the quality and reliability of the Cobolt high performance lasers. The C-FLEX laser combiners are available in three platform sizes (C4, C6, and C8) and can be equipped with up to 4, 6, or 8 Cobolt lasers. C-FLEX can be fully customized, or is available as application-specific configurations. The design allows for full flexibility in the choice of laser technology, ranging from plug and play diode lasers to high power, single frequency diode pumped lasers.

#### The way setups are meant to be

C-FLEX features a common power supply and common interlock (key switch plus remote interlock) for all lasers. It is field upgradeable and customizable. The compact and robust design of the C-FLEX provides excellent long-term stability and outstanding flexibility for your application.

#### **Highest flexibility**

- 375 nm to 1064 nm, up to 1000 mW
- Easy to install and field upgradeable
- High speed modulation, optional AOM (acousto-optical modulators) for fast modulation of DPSS lasers and power regulation
- Fiber coupling with options for multiple outputs
- Optional electromechanical aperture shutters
- User friendly software interface

#### Cant see the wavelength of your choice? Contact our sales team.

#### Applications

- Fluorescence microscopy
- Flow cytometry
- Raman spectroscopy
- Optogenetics
- Photochemistry
- Holography
- Argon-Ion replacement
- Custom Solutions



## VALO Ultrashort femtosecond fiber lasers



#### Typical optical spectrum



#### Typical temporal profile



#### **Fiber lasers**

In general fiber lasers offer a lot of advantages compared to solid state lasers. Due to the waveguide effect of optical fibers the laser systems have very good thermal and vibrational stability and can produce nearly diffraction limited beam profiles. Fiber lasers offer a low cost of ownership and they are essentially maintenance free.

#### Utilizing the nonlinear effects

The VALO Femtosecond Series lasers offer pulse durations of below 50 fs combined with remarkable pulse quality and different energy levels. Due to their unique design, with a perfect interaction of linear and nonlinear effects, an optical bandwidth much larger than standard lasers can be achieved.

The lasers are passively cooled with no fan or water cooling. Their short pulses and resulting high pulse peak powers are especially interesting in applications like multiphoton microscopy since they allow for higher photon conversion efficiency while reducing photodamage and increasing cell viability. The lasers are built into a compact and robust package, making them suitable for integration into industrial screening tools and benchtop analytical instruments.

#### Typical power stability



#### **Applications**

- Multiphoton microscopy
- Optogenetics
- Two-photon polymerization
- Terahertz generation



## **Our laser products –** A wide array of wavelengths

	VA-OSaleson	Unesels Canone	seres cool of the seres	AND ROOS	opologing the state of the stat	A loso seres soleron seres months	olegies stuges capacity capacity	No seven as a seven as	Ne traff	store the second	adde agent						
													0 10 20	50 70	80 10	0 15	0
320 nm	-	-		<u> </u>	-	-		-	~	-	_	320 nm	05-01				
355 nm					•						-	355 nm	05-01		-	or™	
375 nm				•	-	•						375 nm		06-01			
395 nm				•		•					-	395 nm			06-01		
405 nm							•	•			-	405 nm	08-01		Skyra™		
400 mm								-			-	415 nm			06-0	1	
425 nm											-	425 nm			06-0		
425 nm								•			-	445 pm	Sk	wraTM			
450 pm				•		•		•				450 pm				-	C 14/41
450 mm					•		•					450 1111	00.01 04-0	1/06-01		-	C-WA
457 mm				•	•	•	•	•			_	457 1111	08-01 / 5	Skyra™	01 714		C-WA
473 nm				•	•	•	•	•				473 nm		04-01	Skyra™		C-WA
488 nm				•		•	•	•				488 nm	08-01		Skyra™	C-WAV	E / 06-0
491 nm				•	•						_	491 nm		04-01		C-WAVI	= / 05-0
505 nm				•		•						505 nm		06-0	1		C-WAV
515 nm				•	•	•	•	•				515 nm	08-01 / Sk	syra™		04-01 / 06-01	C-V
520 nm				•		•						520 nm		06-0	1		
528 nm				٠		•						528 nm		06-0	01		
532 nm				•	•	•	•	٠	•			532 nm	Tor	™ XS	Skyra™	06-0	)1 / 08-(
553 nm				•		•		•				553 nm	06-01 / Sł	kyra™			
561 nm												561 nm	Sł	(yra™ 0	6-01 / 08-01		04-0
594 nm				•	•							594 nm			04-01		
633 nm				•		•	•	•				633 nm	08-01 SI	kyra™ 06-0	01		
638 nm				•		•	•	•				638 nm	SI	kyra™ 08-0	01		06-01
640 nm					•					•		640 nm					
* Center wavelen	oth selectable b	etween 2-5um	n. tunable 50nm.										0 10 20	50 70	80 10	0 15	.0

\* Center wavelength selectable between 2-5um, tunable 50nm.

\*\* C-WAVE typical powers can reach up to Watt-level in the visible and infrared wavelength range. Compare with typical tuning curve for more details.



## **Our laser products –** A wide array of wavelengths

	VALO Seites	nd beens	asets coloritation	Strike Ball Constant	onthe cool of the strate	allos seies enerchises broute	or serves entrepts coloring theres	November of the second	N <sup>det</sup> co <sup>oft</sup> co	ond lingth of the set								
				ATT.			<b>`</b>					0	10 20	50	70 80	) 10	0 1	150
647 nm	-		~		~				•	~	-	647 nm		Skyra™			06-01	
650 nm		-		-				-			-	650 nm						
660 nm		-			•		•				-	660 nm	Skyra™	1/08-01	04-0	1 / 06-01		
695 nm		-			•		•				-	695 nm	Oryla	01	04-0			
000 mm		-									-	000 mm	-06-					
690 nm		-		•		•					-	690 nm						0
730 nm		-		•							-	730 nm		06-01				
750 nm												750 nm						
760 nm				•								760 nm	06-01					
785 nm				•			•					785 nm						
808 nm				•								808 nm				06-	-01	
830 nm				•			•					830 nm				08-01		
852 nm				•		•						852 nm		06-01				
900 nm												900 nm						
915 nm				•		•						915 nm						
940 nm				•								940 nm						
975 nm				•								975 nm						
1020 - 1080 nm												1020 - 1080 nm					١	/ALO A
1064 nm				•	•		•		•			1064 nm				Tor™ XS		06
1500 nm												1500 nm						
1700 nm*												1700 nm*						
2200 nm*												2200 nm*						
3264 nm*												3264 nm*			Odin™			C-W
3431 nm*												3431 nm*			Odin™			C-W
4330 nm*												4330 nm*			Odin™			
* Center wavele	noth selectable	between 2-5un	tunable 50nm									(	) 10 20	50	70 80	0 10	0 1	150

Center wavelength selectable between 2-5um, tunable 50nm.

\*\* C-WAVE typical powers can reach up to Watt-level in the visible and infrared wavelength range. Compare with typical tuning curve for more details.



## APPLICATIONS Applications

Raman spectroscopy		•	
Life Science:			
- Fluorescence microscopy	•		
- Multiphoton microscopy	•		
- Photoacoustics			
- Flow cytometry			
- DNA sequencing			
Semiconductor inspection			
Quantum research		•	
Interferometry		•	
LIBS			

0

## Laser Solutions for many Applications

We expand the capabilities of your application – both in research as well as industry. As a highly experienced system supplier of technical industrial products, HÜBNER Photonics is developing innovative technology to meet the demands of the global market.



## Applications Lasers for Raman Spectroscopy

The "inelastic scattering of light" was first observed and identified by 
Cobolt – Single frequency and narrow linewidth lasers C.V. Raman in 1928. In 1930, he received the Nobel Prize in Physics C-WAVE - Tunable laser, single frequency for this work. But only in more recent times has it become possible • C-FLEX - Laser combiners to make use of the so-called Raman effect by which the frequency of light changes when it is deflected by molecules. In the past two decades, Raman spectroscopy has developed as a widely applicable method of examination - for fields ranging from material analysis to life sciences applications to point-of-care diagnostics. This has been made possible through the development of compact laser sources, highly sensitive cameras and compact high-resolution spectrometers.

With the lasers from Cobolt, HÜBNER Photonics has one of the widest ranges of products on the market for Raman spectroscopy applications. The Cobolt 08-01 Series of lasers, for example, has been specially developed for Raman applications with integrated Raman filters and optional isolators. Individual lasers can be combined with the C-FLEX laser combiner for added user friendliness. For applications requiring a tunable, single-frequency laser in the visible range, the C-WAVE laser is especially well suited.



## **Applications** Lasers for Life Science

Significant applications in life science requiring lasers include fluorescence microscopy, multiphoton microscopy, photoacoustics, flow cytometry and DNA sequencing. In all of these applications the fluorescence of specific fluorophores or biomarkers is detected, counted • C-FLEX - Laser combiners or imaged, leading to a deeper understanding of biochemistry.

With the lasers produced by Cobolt, HÜBNER Photonics can provide a broad portfolio of high-performance lasers for fluorescence microscopy applications. Not only do these lasers feature a standard integrated clean-up filter, at >70 dB they have one of the best modulation-extinction ratios in the industry. For greater convenience and ease of use, the individual lasers can be combined using the C-FLEX laser combiner. For multiphoton microscopy, the VALO sub 50 fs fiber lasers result in greater efficiency with lower total average power at the sample and therefore allowing longer imaging times - a true advantage.







- VALO < 50 femtosecond pulses</p>
- Cobolt 06-01 Modulated lasers
- Cobolt Skyra<sup>™</sup>– Multi-line laser
- Cobolt Tor<sup>™</sup> Series Nanosecond lasers



Typical pulse profile and beam profile of VALO femtosecond series









### Applications Lasers for Semiconductor Inspection

Semiconductor inspection includes many different techniques, some 

C-WAVE – Tunable laser, single frequency of which use lasers, for the metrology and inspection of semiconduc- • Cobolt - Single frequency and narrow linewidth lasers tor wafers during the manufacturing process.

A few of the inspection techniques used during this inspection process include interferometry, ellipsometry, particle analysis, photoluminescence, fluorescence and Raman spectroscopy. The high performance lasers from HÜBNER Photonics are commonly used for many of these inspection techniques due to their extremely good performance specifications and excellent robustness, meeting the requirements of this demanding application.



### Applications Lasers for Quantum Research

The field of quantum optics deals with the interactions of light and C-WAVE – Tunable laser, single frequency matter. In this still young discipline, photons or "light quanta" - the Cobolt - Modulated diode and single frequency lasers smallest particles of light - are investigated. The aim of this type of research is to gain an understanding of the overall behavior of electromagnetic waves.

The laser sources produced by HÜBNER Photonics play an important role in quantum optics research. This is especially the case with the tunable C-WAVE laser and Cobolt lasers with their high spectral purity and excellent wavelength stability. They provide exactly the flexibility and the precision that are essential for basic research applications.



Exemplary measurement of stabilized wavelength using AbsoluteLambda™.

### Applications Lasers for Interferometric Applications

The term interferometry derives from the word interference. Inter- Cobolt - Single frequency lasers ference is a phenomenon that occurs when two waves of any kind • C-FLEX - Laser combiner come together at the same time and place. Interferometry makes use C-WAVE - Tunable laser, single frequency of interference phenomena for measurement purposes, for example investigation of the flatness of an optical surface.

Even though there are many different types of interferometric based applications, all of them basically operate on the same principle. Two beams are separated and then combined so that they interfere with each other. In order to get interference a highly coherent laser source is required. The more coherent the laser beam, or the longer the coherence length, the finer the detail that can be resolved.

HÜBNER Photonics has one of the widest ranges of single frequency lasers in the industry for interferometric techniques, including holography, Doppler velocimetry and dynamic light scattering.



## **Applications** Lasers for LIBS Laser-Induced Breakdown Spectroscopy

Laser-induced breakdown spectroscopy (LIBS) is an atomic-emission Cobolt Tor<sup>TM</sup> – Nanosecond lasers spectroscopy technique that enables rapid chemical analysis of a wide range of materials ranging from metals, semiconductors, glasses, biological tissues, plastics, soils, thin-paint coating, and electronic materials. LIBS based imaging is extending the capability of this technique even further by providing information on the distribution of critical elements within a sample (as seen below).

The Cobolt Tor<sup>™</sup> Series lasers represent a new class of compact high performance diode-pumped Q-switched lasers that could really boost the trend of extending the use of LIBS systems from laboratory work to industrial applications. A key advantage of the Cobolt Tor™ Series lasers is their unique combination of compact size and performance. The lasers are also manufactured into hermetically sealed packages ensuring robust performance and long lifetime in varying ambient conditions making the lasers highly suitable for integration into demanding industrial applications.









# **Terahertz technology**

Visualizing the invisible

Terahertz (THz) waves penetrate paper, clothing, plastics and many other materials without causing any damage - contact-free and in seconds. This makes THz technology a unique tool for the detection, analysis and hyperspectral imaging of materials. On the basis of the latest research, HÜBNER Photonics is developing highly compact systems for various applications, ranging from mail inspection to industrial quality control. With their plug-and-play features and intuitive user interfaces, HÜBNER Photonics products also deliver outstanding ease of operation. Due to its non-invasive and non-ionizing properties, THz waves are completely harmless for human tissue. Thus, no expensive safety precautions are required for operation.

With the THz imager T-SENSE, small packages and letters can be scanned and hidden contents can be visualized. This makes it possible to detect dangerous materials such as drugs or explosives before opening the mail.

The Prism Award winning THz spectrometer T-COGNITION has been developed in cooperation with the Fraunhofer Institute for Physical Measurement Techniques (IPM) and is an advanced system for the automated identification of hidden drugs or explosives inside letters or small parcels within seconds on the basis of their characteristic spectroscopic fingerprints. The systems help to improve security in mailrooms, at customs stations or at official offices - without any health risk for the users.

The T-SPECTRALYZER product family comprises turn-key THz spectrometers for scientific research as well as industrial non-destructive testing (NDT) applications. The variety of available options provides highest flexibility to make the T-SPECTRALYZER fit to the customers' individual applications.



T-SENSE user interface with different filters side by side. All parameters are intuitively customizable.



## THz waves – Combining characteristics

THz waves are a part of the electromagnetic spectrum in the frequency range of 0.1 THz to 10 THz. Many non-conductive materials such as plastics or PVC, ceramics, paper or clothing are nearly transparent in the THz frequency range, a characteristic feature known from GHz waves. At the same time, THz waves are absorbed at varying degrees by substances such as drugs, explosives and many other materials, providing spectroscopic information about the sample material, similar to a spectrum produced using optical wavelengths. These absorption characteristics can be used to produce a sort of spectroscopic fingerprint, making it possible to identify various substances, even when they are concealed in a letter or parcel, for example.



## Applications

#### **Public security**

- Mail inspection by hyperspectral imaging for detection of hazardous goods
- Identification of drugs and explosives using Time-Domain-Spectroscopy (TDS)

#### Non-destructive testing (NDT)

- Detection of defects and cavities in components (hyperspectral imaging)
- Identification of substances through plastic pipes or tubing and other types of packaging using Time-Domain-Spectroscopy (TDS)
- Determination of the thickness of layers of multi-layered systems for identification purposes
- Analysis of substances in powder and tablet form
- Analysis of gases
- Investigation of moisture distribution
- Distinction between crystal and amorphous structures







## **T-SENSE**

#### The THz imager for safe and secure mail inspection

Mailrooms in companies and in public agencies are particularly vulnerable to attacks via postal deliveries. Often the personnel at these facilities are not protected against the dangers of such attacks. To reduce these risks, HÜBNER Photonics has developed an innovative table-top device for the inspection of letters and small packages. The THz imager T-SENSE visualizes objects and hazardous substances that may be concealed in postal items.

Unlike x-rays, there is no health risk with the THz imager T-SENSE. Typical safety precautions for conventional imaging devices are not needed. And thanks to its intuitive graphic user interface, the systems are very simple to operate. T-SENSE facilitates effective and safe postal inspection, helping to keep people safe and secure.

## **T-COGNITION**

#### Identifying hazardous substances

The THz spectrometer T-COGNITION is an extremely effective security technology. With reliability and precision, it identifies hidden drugs and explosives in letters and small packages without the necessity of handling or opening the item in question. Within seconds, T-COGNITION identifies the spectroscopic fingerprint of the hazardous substance or material by comparing the data with its own database. This system enhances work safety in prisons, at custom controls, at authorities, in companies and embassies, to name but a few.

## **T-SPECTRALYZER**

#### The all-round spectrometer

The T-SPECTRALYZER products have been developed for quick use applications and for routine measurements in daily analytical work. The product family comprises turn-key systems operating either in transmission or reflection geometry or both simultaneously as well as fiber-coupled solutions, which provide full flexibility to adapt the spectrometer to the customer's application. The systems only require a normal mains connection and are ready for use without any additional infrastructure. Thanks to the most modern technology, the THz spectrometers function without any additional cooling or external gas supply. This allows for quiet and very economical operation. Individual add-on modules and the intuitive user interface support the recording, processing and export of measurement results. Not least, the automatic user-configurable report functionality makes T-SPECTRALYZER the most user-friendly THz spectrometer.



#### **HÜBNER** Photonics

#### HÜBNER Photonics GmbH (Sales in Germany, Switzerland and Austria)

Wilhelmine-Reichard Strasse 6 34123 Kassel, Deutschland

Phone: +49 561 994 060 - 0 Sales: +49 561 994 060 - 11 E-mail: info.de@hubner-photonics.com

#### Cobolt AB (Sales in Norway, Sweden, Finland and Denmark)

Vretenvägen 13, SE-171 54 Solna, Sweden

 Phone:
 +46 8 545 912 30

 Fax:
 +46 8 545 912 31

 E-mail:
 info.se@hubner-photonics.com

## VALO Innovations, a part of HÜBNER Photonics

Hollerithallee 17 30419 Hannover, Deutschland

Phone: +49 511 260 390 70 E-mail: sales@hubner-photonics.com

#### HÜBNER Photonics UK (Sales in UK & Ireland)

1a Northcote Avenue Isleworth, Twickenham Middlesex, TW7 7JH, London United Kingdom

Phone: +44 7359440871 E-mail: info.uk@hubner-photonics.com

#### **Direct sales offices**

HUBNER Photonics Inc. (Sales in USA, Canada & Mexico)

2635 North First Street, Suite 202, San José, California, 95134, USA

 Phone:
 +1(408)708 4351

 Fax:
 +1(408)490 2774

 E-mail:
 info.usa@hubner-photonics.com

#### www.hubner-photonics.com

#### **Sales representatives**

Australia, New Zealand Warsash Scientific Pty Ltd www.warsash.com.au

Benelux Laser 2000 Benelux CV www.laser2000.nl

Brazil Photonics Instrumentos www.photonics.com.br

China DynaSense Photonics Co. Ltd. www.dyna-sense.com

Estonia, Latvia and Lithuania Optek Ltd www.optek.lv

France Optoprim www.optoprim.com

India Spectral Instrument System www.spectral-instruments.com

Israel Lahat Technologies Ltd www.lahat.co.il

#### Italy

Crisel Instruments www.crisel-instruments.it

Japan Kantum Electronics Co Ltd www.kantum.co.jp

Pneum Co, Ltd Japan www.pneum.co.jp

Japan Laser Corporation www.japanlaser.co.jp

Poland Amecam www.amecam.pl

Singapore, Malaysia and Thailand Wavelength Opto-Electronic (S) Pte Ltd www.wavelength-tech.com

South Korea BM Laser Solutions www.bmlaser.co.kr

Spain and Portugal Laser Technology SI www.laser-technology.com

Taiwan Tayhwa Technology Co Ltd www.tayhwa.com.tw

## www.hubner-photonics.com







www.hubner-photonics.com