High-resolution High-Speed LIBS Imaging

A.Tercier^{1,2}, E. Vasileva^{3*}, C. Alvarez-Llamas¹, C. Fabre⁴, S. Hermelin¹, F. Trichard², C. Dujardin¹, V. Motto-Ros^{1,2**}

Corresponding authors: * Elena Vasileva, Ph.D, Product Manager at Hübner Photonics.com, ** Vincent Motto-Ros, Associate Professor, Insitut Lumière Matière UCBL vincent.motto-ros@univ-lyon1.fr

Motivation

Laser Induced Breakdown Spectroscopy (LIBS) is becoming increasingly prevalent as an analytical technique for material content due to its high sensitivity and the capability to detect any element of the periodical table with nearly no requirements or limitations on sample dimensions, consistency, and/or surface quality. LIBS imaging has expanded the capabilities of LIBS technique even further [1]. LIBS imaging enables the reconstruction of elemental or mineralogical maps which show the distribution of the elements or minerals within an analyzed sample. Typically, excitation sources, such as lasers, with a repetition rate between 10-100 Hz are used for LIBS analysis. These systems have traditionally provided a sufficient balance of sampling frequency and laser pulse energy. However, for high-resolution µmscale LIBS images of large sample areas (several cm²) a higher repetition rate is required to avoid unreasonably long acquisition times. This work demonstrates an approach to reduce the acquisition time for high resolution µ-LIBS imaging by using a laser operating in the kHz frequency range [2]. The faster laser repetition rate enabled the development of a **µ-LIBS imaging microscope** capable of resolving images with about **10 µm resolution at < 20 min/cm²**. For the first time, **images of such high** resolution showing detailed elemental distribution within an analyzed sample have been achievable. The capability of the system to perform fast LIBS analysis over large samples and revealing the spatial distribution of elements within the analyzed area opens opportunities for application within a wide variety of fields of research including **biomedical and geological material analysis**, as well as **in industrial** settings such as mining.

Cobolt Tor[™] XE


Compact air-cooled ns-pulsed laser



Cobolt Tor[™] XE


Pulse energy (µJ) @ 1064 nm	500 ± 50
Pulse length (ns)	1 - 3
Spatial mode (TEM00), M ²	< 1.3
Pulse jitter (μs)	< 2
Power stability (8 hrs ± 3°C)	< 3
Beam diameter (mm)	0.7 ± 0.1
Repetition rate (kHz)	up to 1 k
Trigger source	External,


In LIBS systems, the choice of the laser used for material ablation plays an important role. LIBS measurement repeatability can be linked to laser performance as the plasma generation relies on nonlinear processes. Any fluctuations of laser parameters such as **pulse energy** (or **pk-to-pk stability**), pulse duration, and beam quality, could produce fluctuations in ablated mass, plasma temperature, and electron density, which would be reflected in the emission spectrum and could lead to analytical inaccuracies. Laser performance is even more critical in LIBS imaging, where raster scanning of the sample typically means using a single laser pulse per acquisition point. The quality of the beam M² is associated with the ability to focus the beam into a spot which defines the spatial resolution of the system, thus, the min.pixel size of an image. The repetition rate of the laser defines the maximum acquisition/scanning speed available. The footprint of the laser and the ability to air-cool it, can be big advantages for system integration.



The proposed µ-LIBS imaging system is a custom-made LIBS µ-microscope capable of resolving images with about 10 μ m resolution at < 20 min/cm². The system is based on a compact nanosecond Nd:YAG Qswitched 1064 nm laser (Cobolt Tor[™] XE), with a pulse energy of 0.5 mJ and a pulse duration of 2.5 ns. The laser beam is focused down to 10 µm spot. The plasma emission lines are recorded with a spectrometer (Kymera 328i from Andor Technology) coupled to a sCMOS camera (iStar from Andor Technology).

High-resolution 10 µm high-speed < 20 min/cm² LIBS image





a) Photo of a sample placed under the µ-LIBS microscope.

The high-speed, high-resolution µ-LIBS imaging results above were enabled by gated sCMOS technology and an ultra-stable pulse kHz laser. To demonstrate imaging capabilities of the µ-LIBS imaging microscope over large sample areas, a rock sample with an area of about 6 cm² was used. The sample was polished with SiC paper under water solvent to obtain a planar and clean section. The image shows the distribution of silicon (Si) for quartz, lead (Pb) for galena, calcium (Ca) and magnesium (Mn) for carbonate, and iron (Fe) for pyrite in a sample with 10 µm resolution. The reconstructed µ-LIBS image was completed in less than 2 hours (< 20 min/cm²), to put it into perspective, similar image would take 20 hours to complete with a laser operating at 100 Hz [3].

kHz , Internal, Gated



b) Reconstructed µ-LIBS image showing the distribution of silicon (Si), lead (Pb), calcium (Ca), magnesium (Mn), and iron (Fe) in a sample.

iStar from Andor Technology Fast gated camera based on sCMOS sensor technology iStar sCMOS (13 µm pix [2x2 bin]) **Pixel read noise** 5.2 e-(full frame max. rate) Gain to overcome ~X31 sensor noise floor Pixel well depth 120,000 e-Dynamic range ~3,870:1 APS with column A/D 전전전전 왜 오래 오래 오래 오래 그 A/D A/D A/D A/D A/D READOUT Horizontal Scan Circuit vIOS allows for low readout noise at high readout rates allowing for fast data acquisition while maintaining high dynamic range

Latest generation fast gated cameras based on sCMOS sensor technology provide significant enhancement in terms of high acquisition rates and simultaneously high dynamic range compared to CCD, Interline or EMCCD-based gated detectors.


Summary

- equates to $< 20 \text{ min/cm}^2$.
- can specifically.

References

HÜBNER Photonics





A high-speed high-resolution µ-LIBS imaging system, µ-LIBS microscope, capable of reconstructing elemental compositional images with 10 µm resolution over a sample area of 6 cm² in a dramatically reduced time of < 2 h which

A drastically reduced time for material content analysis of relatively large sample area with µmprovide unique capabilities and opens new perspectives for LIBS imaging and LIBS analysis in many application fields, and in geological material analysis

Acknowledgement

This work was partially supported by the French region Rhônes Alpes Auvergne (Optolyse, CPER2016), the French "Agence Nationale de la Recherche" (ANR-22-xxx "MEMOar" and ANR-20-CE17-0021 "dIAg-EM), and the French government "Plan de Relance" and "France 2030" (ANR-22-PEXD-0014, "Libelul") under the DIADEM program managed by the "Agence Nationale de la Recherche" (ANR-22-PEXD-0014, "Libelul"). In addition, we gratefully acknowledge Jean Michel Laurent, Shayne Harrel, Antoine Varagnat from Andor Technologies.

[1] Gardette V, Motto-Ros V, Alvarez-Llamas C, Sancey L, Duponchel L, Busser B, Analytical Chemistry 2023 95 (1), 49-69, DOI: 10.1021/acs.analchem. 2c04910

[2] V. Gardette V, Motto-Ros V, Alvarez-Llamas C, Sancey L, Duponchel L, Busser B. Laser-Induced Breakdown Spectroscopy Imaging for Material and Biomedical Applications: Recent Advances and Future Perspectives. Anal Chem. 2023 Jan 10;95(1):49-69. DOI: 10.1021/acs.analchem.2c04910. PMID: 36625118.

[3] Jolivet L, Leprince M, Moncayo S, Sorbier L, Lienemann C-P, Motto-Ros V, Review of the recent advances and applications of LIBS-based imaging, Spectrochimica Acta Part B: Atomic Spectroscopy. 151 (2019) 41-53. DOI: 10.1016/j.sab.2018.11.008.